

THE MINI GIANT



iBon Mini &
iBon Mini-SE

iBon Mini

Brushless Mini Endo Motor

➤ Smaller Contra-angle

48mm x 15mm x 15mm

Allows unobstructed view of operating field

➤ Brushless Motor



Wider speed range: **150-1500rpm**
Lower noise




Stable torque
control

➤ Light and Compact

- **99g/190mm**

- Lightweight

- No strain on operator even after prolonged use



Type	Brushless & Cordless
Speed/Torque	1500 rpm Max/4.0 Ncm Max
Weight/Length	99g/190mm
Contra-Angle	6:1
Contra-Angle Size	48mm*15mm*15mm
Rotation	Clockwise/Anti-Clockwise/Reciprocating
Reciprocation	CCW/CW Angles Adjustable
Memory	10 Settings
Charging	Type-C




Integrated
Bluetooth

MINI BUT MIGHTY




↘ Light and Compact

- 99g/190mm
- Lightweight
- No strain on operator even after prolonged use




↘ Slimmer Contra-angle

- 48mm x 15mm x 15mm
- Allows unobstructed view of operating field



iBon Mini-SE

Cordless Mini Endo Motor



Integrated
Bluetooth

Type	Brush & Cordless
Speed/Torque	800 rpm Max/4.0 Ncm Max
Weight/Length	99g/190mm
Contra-Angle	1:1
Contra-Angle Size	48mm*15mm*15mm
Rotation	Clockwise/Anti-Clockwise/Reciprocating
Reciprocation	CCW/CW Angles Adjustable
Memory	10 Settings
Charging	Type-C

MINI BUT MIGHTY



Distributed By.



Bondent Group
Email: international@bondent.com



@Bondent Group